MD

FOSHAN CICADA DENTAL INSTRUMENT CO., LTD
B5-2F, Guangdong New Light Source Industrial Base,
Nanhai District, Foshan, Guangdong, 528226, P.R. CHINA

EC REP

Lotus NL B.V.

Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague – NETHERLANDS

CADENCE

2 bis, chemin du loup
93290 Tremblay-en-France - FRANCE

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1. Recommendation of Use

Before using, please read this instruction manual carefully and then keep it in a safe place.

1.1 Intended Use /Indication for use

LED Curing Light is used for light curing polymerization of dental composites, luting materials, cements and other light cured materials.

1.2 Contraindication

Show caution when using this equipment with patients with heart disease, pregnant women or children.

1.3 Warning

- If the curing light works for 40s continuously, the temperature of the top of optical fiber may reach ~56 °C.
- Do not modify this equipment without authorization of the manufacturer.
- The adapter should be connected to a socket which is easy for operator to reach.
- After a working cycle, the operator can press the ON/OFF switch to start another working cycle. Stop operating if the equipment starts to get noticeably hot, do not restart until the equipment has cooled down. We recommend not to perform more than 10 working cycles immediately after the other.
- The light radiation produced by this type of unit is dangerous and must not be pointed at eyes. The light produced by this equipment must be directed only at the zone to be treated in the oral cavity.
- Do not use the unit near a heat source. Do not use solvents, detergents or flammable products, which may damage the unit or cause a short-circuit.
- The device is not waterproof.
- Although the device conforms to standard of TEC 60601-1-2, it may still cause or be affected by EMI of other devices. If in case it happens, separate the interfering devices from each other, re-orient cabling, or plug device into a different power outlet. If it does not correct the problem, contact the manufacturer for technical support.
- Do not use the system in close proximity to strong magnetic fields (such as found near an MRI).
- The device is not provided sterile and is not intended to be sterilized when used. To prevent cross contamination, users must cover the device with a FDA cleared barrier sheath sleeve. It should be changed between each patient.


1.4 Precautions

- **Charge the battery for at least 4 hours before first use.**
- The top of the unit can be rotated through 360 degrees, it must not be removed when in operation.
- During operation, the blue light should be aimed straightly at the composite resin to ensure the effect of solidification.
- Avoid aiming the blue light at eyes directly.
- Please use the power adapter and lithium battery which is designed and supplied by our company. Using power an adapter and a lithium battery which is designed and supplied by other manufacturers may cause potential dangers to the operator and the patient.

- It is forbidden to use metal or other conductors to touch the main unit and the charging point of pedestal because it may burn the internal circuit or make the lithium battery short circuit.

2. Structure and components

LED-3200 curing light is composed mainly of a high power LED and a main unit.

P/T button

Press the P/T button for about 2 seconds to select the operation mode, the device provides four modes, Normal Mode, Pulse Mode, Strong Mode, and Super Mode.

- Normal Mode (P1): Mode with 5, 10, 15, 20 seconds exposure time.
- Pulse Mode (P2): Mode with 5, 10, 15, 20 seconds exposure time option.
- Strong Mode (P3): Mode with 3, 6, 9 seconds exposure time option.
- Super Mode (P4): Mode with 1, 2, 3 seconds exposure time option.

ON/OFF button

Press the ON/OFF button to start the curing light according to the set modes, and the main unit will produce a “Bi” sound. It can also stop the curing light anytime.

Display screen

When the device is turned on, the work mode will be shown on the display screen:

Normal mode (shown as 'P1'); Pulse mode (shown as 'P2'); Strong mode (shown as 'P3'); Super mode (shown as 'P4'); Press Mode button to switch the model for P1; P2; P3; P4 and after select the mode, press Time button to select the exposure time.

3. Technical specifications

Product Name	LED Curing Light
Model	LED-3200
Intended user	Dentist
Intended use site	Hospital or dental office
Light source	High-power blue light LED
Wave length	385 nm to 515 nm

Light intensity	800 mW/cm ² to 3200 mW/cm ²
Depth of solidification	> 2 mm
Work Mode	<ul style="list-style-type: none"> - Normal Mode (P1): 5, 10, 15, 20 seconds, 800-1000 mW/cm² - Pulse Mode (P2): 5, 10, 15, 20 seconds, 800-1000 mW/cm² - Strong Mode (P3): 3, 6, 9 seconds, 1200-1400 mW/cm² - Super Mode (P4): 1, 2, 3 seconds, 3000-3200 mW/cm²
Power source	Rechargeable Lithium battery, 3,7 V / 2000 mAh
Charger	Power Input: 100 V to 240 V, 50 Hz/60 Hz Input Power: 10 VA Power Output: DC 5V / 1A
Operation Mode	Intermittent operation
Protection type against electrical shock	Class II
Protection degree against electrical shock	Type B 
Protection against harmful ingress of water or particular matter ordinary equipment	IPX0
Safety in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide	Not suitable under this condition
Weight	136 g
Dimension	24 (W) x 240 (H) mm

4. Installation

- The top of the unit can be turned 360 degrees to both left and right, it is forbidden to remove when working. Insert the top of the unit into the main unit in rotation mode.
- When the battery needs to be charged, connect the plug of the adapter into the ac 100V- 240V power supply. Then connect the output plug of the adapter to the de DC 5V input jacket of the charger, then put the main unit into the charger and connect it to the charger station terminal.

5. Operation

- The device is not provided sterile and is not intended to be sterilized when used. To prevent cross contamination, users must cover the device with a FDA cleared barrier sheath sleeve before use. It should be changed between each patient.
- Press the P/T button, following four modes are available: Normal Mode: screen show P1. Press the P/T button to select of time which could be 5, 10, 15, 20 seconds.
- Pulse Mode: screen show P2. Press the P/T button to select of time which could be 5, 10, 15, 20 seconds.
- Strong Mode: screen shows P3. Press the P/T button to select of time which could be 3, 6, 9 seconds.
- Super Mode: screen shows P4. Press the P/T button to select of time which could be 1, 2, 3 seconds.
- During the operation, aim the top of the main unit, at the correct position, press the ON/

OFF button and the main unit will produce a “Bi” sound, the curing light radiates blue light and starts working according to the set modes meanwhile. It starts counting down and will produce a “Bi”, it stops working when counting down to 0. At the end of the work, the digital display returns to the set work mode.

- During operation, the blue light can be stopped by press the power button at any time.
- After a working cycle, the operator can press the ON/OFF switch to start another working cycle. Stop operating if the equipment starts to get noticeably hot, do not restart until the equipment has cooled down. We recommend not to perform more than 10 working cycles immediately after the other.
- Low power detective circuit is fixed inside of the main unit when low power is detected, the display screen of main unit will wink, please charge in time.
- When the battery needs to be charged, connect the plug of the adapter into the ac 100 V ~ 240 V power supply. Then connect the output plug of the adapter to the DC 5.0 V input plug of the charger station. Put the main unit to the charging point of the charger station, the display screen turned to scrolling display, and the curing lights starts charging.
- This equipment will turn off automatically if there is no any action within 7 minutes; You can turn it on by press any button.
- The depth of solidification of composite is no less than 4mm per 10 second.

6. Clean and Disinfection

In order to prevent cross-infection, it's commended to do following process before and after each operation.

1. Turn off the device and take the main unit from the charger.
2. Take off the cover sleeve after operation.
3. Wipe the main unit with a soft cloth moistened with purified water.
4. Let the device air dry.

7. Daily maintenance

- This equipment does not include the self-maintainable spare parts.
- The maintenance of this equipment should be taken by the appointed professional or special repair shop.
- Users can change the lithium battery. Please use accessories which are designed and supplied by our company, contact with the local dealer or our company if you want to buy. Accessories designed and supplied by other manufacturers may cause potential hazards to the unit or other damage.
- The accessories of the product should be cleaned with clean water or neutral sterilized liquid. Do not soak. Do not use highly volatile and diffluent solvent to clean this equipment, which can cause the signs on the control panel to fade.
- Please clean the resin remained on the top of the main unit after using to avoid infecting the life-spanor solidified effect.
- The device cannot be maintained during operation, and it is suggested to maintain it once a month.

8. Trouble Shooting

Fault	Possible cause	Solutions
No indication, no response.	<ol style="list-style-type: none">1. Battery is out of power.2. Faulty of battery.3. The main unit battery protection system works.	<ol style="list-style-type: none">1. Charge the device.2. Send to after service for repair.3. Place the main unit into the socket on the charger for activation.
“Er” shown on the screen.	<ol style="list-style-type: none">1. System error.2. Faulty of main unit.	Send to after service for repair.
Wink shown on the screen.	Low battery.	Reconnect the charger, if “Lo” show again after 15 minutes please send to after service to repair.
Light intensity is weak.	There is resin on top of the main unit.	Remove the resin.
The device is not charging when the adapter is connected.	<ol style="list-style-type: none">1. The adapter is not well connected.2. Faulty of adapter or incompatible.3. The charging unit is dirty.	<ol style="list-style-type: none">1. Reconnect.2. Change the adapter.3. Clean with alcohol.
The mode indicator twinkles when charging.	<ol style="list-style-type: none">1. Low voltage.2. Short circuit of the battery.	<ol style="list-style-type: none">1. Back to normal after 15 minutes charging.2. Send to after service for repair.

If such handlings are completed, the device still cannot work normally, please contact the special maintenance shop or our company.

9. After service

From the date this equipment has been sold, based on the warranty card we will repair this equipment free of charge if it has quality problems, please refer to the warranty card for the warranty period.

10. Storage and transportation

- This equipment should be handled carefully, kept away from shaking point, installed or stored at shadowy, dry, cool and ventilated places.
- Don't store it together with articles that are combustible, poisonous, caustic and explosive.
- This equipment should be stored in the environment where the relative humidity is 10%-93%, the atmosphere pressure is 70 kPa to 106 kPa and the temperature is -20° C to +55 °C.
- Excess impact or shake should be avoided during transportation.
- Don't mix it with dangerous articles during transportation.
- Keep it away from sun or snow or rain during transportation.

11. Environmental protection

Please dispose according to the local laws.

12. Symbol descriptions

The following symbols may appear in this manual, on the LED Curing Light LED-3200, or on its accessories. Some of the symbols represent standards and compliances associated with the LED Curing Light LED-3200 and its use.

Date of manufacture

Manufacturer

Serial number

Type B applied part

Follow instructions for use

Direct current

DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary

The degree of avoid ingress of water or particulate matter into ME equipment

13. EMC Declaration

Guidance and manufacturer's declaration -electromagnetic immunity

The LED Curing Light LED-3200 is intended for use in the electromagnetic environment specified below. The customer or the user of the LED Curing Light LED-3200 should ensure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 KV contact ± 18 KV air	±6 KV contact ± 18 KV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient /burst IEC 1000-4-4	±2 KV for power supply lines ± 1 KV for input /output lines	±2 KV for power supply lines ± 1 KV for input /output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 KV differential mode ±2 KV common mode	±1 KV differential mode ±2 KV common mode	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply inputlines IEC 61000-4-11	< 5 % UT (>95 % dip in UT) for 0,5 cycle	< 5 % UT (>95 % dip in UT) for 0,5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the LED Curing Light LED-3200 requires continued operation during power mains interruptions, it is recommended be powered from an uninterruptible power supply or a battery.
	40 % UT (60 % dip in UT) for 5 cycles	40 % UT (60 % dip in UT) for 5 cycles	
	70 % UT (30 % dip in UT) for 25 cycles	70 % UT (30 % dip in UT) for 25 cycles	
	< 5 % UT (>95 % dip in UT) for 5 sec	< 5 % UT (>95 % dip in UT) for 5 sec	
Power frequency (50/60 Hz)magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration -electromagnetic immunity

The LED Curing Light LED-3200 is intended for use in the electromagnetic environment specified below. The customer or the user of the LED Curing Light LED-3200 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test I/cv I	Compliance level	Electro magnetic environment-guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the LED Curing Light LED-3200, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter

Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>Recommended separation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz to 800MHz $d=2.3\sqrt{P}$ 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> 
-----------------------------	-------------------------	-------	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF, transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the LED Curing Light LED-3200 is used exceeds the applicable RF compliance level above, the Medical LED Curing Light LED-3200 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the LED Curing Light LED-3200.
 Over the frequency range 150 KHz to 80 MHz, field strengths should be less than [v1] V/m.

Guidance and manufacturer's declaration -electromagnetic emissions

The LED Curing Light LED-3200 is intended for use in the electromagnetic environment specified below. The customer or the user of the LED Curing Light LED-3200 should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The LED Curing Light LED-3200 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The LED Curing Light LED-3200 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communications equipment and the Medical LED Curing Light LED-3200

The LED Curing Light LED-3200 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Medical LED Curing Light LED-3200 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the LED Curing Light LED-3200 as recommended below, according to the maximum output power of the communications

Rated maximum outputpower of transmitter w	Separation distance according to frequency of transmitter m		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters(m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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