Shanghai International Holding Corp. GmbH
Eiffestrasse 80, 20537 Hamburg, GERMANY

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Remarque importante

Autoclave

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Attention, consulter les documents joints

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Certifié MMD93/42/EEC

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Fabricant

TURBINE TR-740L

User's guide

Please read this user's guide carefully and keep for future reference.

Description

Push button, 4 point water spray, double air tracks, quick coupling light guide vitreous, illumination ≥ 25000 LUX. Voltage is $3.3V \pm 0.1V$, longevity ≥ 70 hours.

Good performance; high durability; can accomplish multiple functions such as dental implant and therapy; acts as necessary equipment for general dental treatment and for dentist's daily working.

⚠ Attention!

- For safe use, follow all instructions carefully.
- For use by dental professionals only.
- Unauthorised dismantling, reconstruction or modification is strictly forbidden.
- Please run point testing outside oral cavity beforehand. If you feel any anomalies, stop using immediately and contact the after-sales centre.
- Using a vibrating or bent bur, or a hand-piece with a worn bearing, will cause a high noise. Continued high noise will harm hearing. Do not use a vibrating or bent bur. If high noise appears during use, please replace the shaft core.

- Do not use a bent, cracked, deformed or unofficial bur. By doing so, it may snap or fly-off which could cause injury or accident.
- Please ensure the bur chuck is clean. If there is any contamination inside, it may cause bur clamping to weaken or shaft core vibration.
- Please run the bur using the rotation speed as specified by the bur manufacturers.
- Do not drop to avoid breakage.
- Please load a bur or a standard stick securely before using.
- Please wait for rotation to stop completely before loading/unloading.

Working conditions

- Air pressure 0.25-0.30MPa, must be filtrated.
- Use high-speed dental bur with shaft diameter 1.59-1.60mm and roughness +/- 0.8 μm .
- The amount of hand-piece spray or water spray should be able to be adjusted.
- Hard glass fibre can be repeatedly sterilised in high temperature and pressure.

Product properties

- Bur holding strength 22-45N.
- Noise 55dB.
- At 25-30 atmospheric pressure: standard hand-piece more than 300,000RPM, big torque hand-piece more than 280,000RPM.

Components

Consists of head, hand-piece cover, 6-hole quick coupling (6-hole fibre optic illuminated hand-piece) and cartridge.

Environment-friendly fibre optic hand-piece tube

If the dental unit doesn't have the optical fibre system, please insert our optical fibre tubing and electricity board in the dental unit. It can then use our optical fibre hand-piece, and include the high speed hand-piece which is compliant with international standards.

Fig. 1

Electricity board

- Power in: 12-24V
- Power out: 2.5-3.5V (adjustable)
- Output delay: 0-12s (adjustable)
- Linkage air pressure: $P > 0.2\text{MPa}$
- Control ways: 2 ways

Installation

- Connect power in line with POWER IN 12-24V end.
- Electricity lines from the tube connect with the power out end (L1 or L2) of the electricity board.
- Connect the electricity board and the air in tube with a T-Connector.
- Air1 corresponding to L1, and Air2 to L2.
- Connect the electricity; step on the foot controller; the indicator light is on; fibre optic illuminated hand-piece works. See figure 3 for more detail.

18

Fig. 2


Fig. 3

Notice:

- A: work air (air in)
- W: water in
- D: air return
- C: pulverisation

Replacement of bulb (optical hand-piece) (6-hole) Fig. 4

- Unplug the 6-hole quick coupling.
- Unscrew the tip cap, remove the burned out bulb.
- Replace with a new one, making sure the pins have good contact with the metal strips.
- In case of bad contact, raise the pins (with small sharp tool) before inserting bulb.
- Use the wrench to remove the old LED bulb and change the new LED bulb.

Connect the quick coupling

- Insert quick coupling
- Insert quick coupling into socket connection of the high speed hand-piece, pull it backwards to check it's in tightly and firmly.
- Extract the quick coupling
- Hold the quick coupling by the decorative pattern, to extract.

Notice: the quick coupling must be fixed tightly and be connected correctly before use.


Fig. 4


Fig. 5

19

Bearing detachment Fig. 6.1

If the ball bearing is damaged but its outer ring hasn't dropped off, put the "bearing outer ring detach tool" between the bearing and the rotor. Put the screw bolt outside, then rotate the screw bolt of tool clockwise until the inner ring is removed. If the ball bearing is damaged and its outer ring has dropped off, insert the notch of the "bearing outer ring detach tool" at the inner ring path of ball bearing, then rotate the screw bolt of tool clockwise until the inner ring is removed.


Fig. 6.1

Bearing installation Fig. 6.2

When pre-installing the bearing, please insert the bearing on the shaft first, with the dustproof shield facing outside, then put it on the bearing installation tool. Knock at the upper side of inner ring of bearing gently with the push axle, ensure the bearing is placed in the middle. Then press the head of push shaft in order to make the bearing is closed to the turbine.

Note: the "bearing outer ring detach tool" and "bearing installation tool" are not supplied. Please contact supplier if the tools are required. Equipment can only be repaired by the manufacturer or authorised agent.

Fig. 6.2

Application

For dental operation (drilling, polishing) with high speed dental bur.

Usage

Replacement of cartridge Fig. 7

- Take out the old cartridge
- Insert the standard test bar in the hole of cartridge chuck, then turn the cover remover anti-clockwise and remove the cartridge cap.
- Push the standard test bar to take out the cartridge.

Installing cartridge

- Open the cartridge cap, then take out the old cartridge and insert the new cartridge.

Note: when replacing head cap, don't use the cover remover first because the threads are pinpoint and can be easily damaged.

Fig. 7

Proper bur

- Do not use bent, damaged, or non-concentric burs. Use burs of ISO-standard \varnothing 1.59-1.60mm.

Change of bur Fig. 8

- Press-cards type (push button bur-exchange type)
- When unloading the bur, press the centre of the head hard with your thumb, whilst simultaneously taking out the bur with the other hand, as shown in picture
- When loading the bur, insert it into importing parts of the chuck hole (about 1-2mm). Then press the centre of the head hard with your thumb, whilst simultaneously pushing the bur into the end with your other hand.

Fig. 8

Replacing the O-rings (quick coupling) Fig. 9.1

Worn or damaged O-rings may cause air or water to leak.

When water is found in the exhaust air or water leak from the coupled point, replace four O-rings as shown.

Replacement O-ring Set


Fig. 9.1

Non-Retraction valve (quick coupling) Fig. 9.2

To prevent oral fluids being drawn into the water line, quick coupling has Non-Retraction valve. If water is beginning to leak from hand-piece, replace Non-Retraction valve.

- Remove the back-end gasket.
- Replace Non-Retraction valve.


Fig. 9.2

Clean-head system

- Scrub dirt and debris from the hand-piece.
- Rotate the hand-piece and immerse half of hand-piece head. Rotate and intermittently stop the hand-piece for about 2 to 3 seconds (about 3 times).
- Stop the hand-piece and wipe it dry. If the dirt couldn't be removed from the hole, remove the cartridge and clean by brush.

Cleaning - bulb Fig. 10

Wipe clean the cellular optical bulb entry point with an alcohol-immersed cotton swab. Remove all debris and oil.

⚠ Notice!

Use only clean water to maintain the clean-head system.

Fig. 10

Maintenance

- Daily lubrication to the bearing is essential for extending its using life.
- After daily use, spray lubricant for 1 to 2 seconds by inserting lubricant nozzle all the way into drive air tube.
- If workload is heavy, lubrication should take place half way through the day and after daily use.
- Lubricate before autoclave.

Autoclave procedure

After every use, autoclave procedure should be taken as follows:

- Remove bur from hand-piece and clean hand-piece by scrubbing with a toothbrush and then use alcohol swab ball to clean hand-piece.
- Lubricate hand-piece.
- Place the hand-piece in an autoclaving pouch.
- Sterilize for 18 minutes at 134°C.

Storage and transportation

Environmental temperature: -40°C to +70°C. Relative humidity 10% to 80%.
The scope of air pressure: 50 kPa to 106 kPa.

Fig. 12

Troubleshooting

Trouble	Possible cause	Solution
Loud noise, low rotation speed, cutting force decrease or hand-piece failing to run	Ball bearing damages	Replace ball bearing
Hand-piece failing to spray mist	Spray hole blockage	Clean with probe
Hand-piece water leakage	O-ring or sealing washer aged	Replace aged parts
Normal noise but low rotation speed	Low air pressure	Adjust air pressure
Bur drop-off or fail to engage bur	Non-standard bur or chucking system damages	Replace new bur or send it to maintenance centre
Bur wobbling, low cutting force	O-ring or ball bearing damages	Replace spare parts

Our service centre can offer technical assistance to you.

Guarantee

Our company grants the user a 12 months guarantee for its complete product range, except ball bearings, from the date of invoice issued. Maintenance over the term of guarantee will be at the customers charge.

Our company will not be responsible for damage or injury resulting from:

- Excessive use.
- Improper manipulation of the product, or modification to the product carried out by persons not authorised by our company.
- Failure to follow instructions to install, operate and maintain the hand-piece.
- Damage of chemical, electrical or electrolysis due to improper autoclaving and/or storing.
- Improper working pressure.

Standard symbols

On the instrument label appear standard symbols as follows:

Important notice

Autoclave

Class 2 equipment

Attention, consult accompanying documents

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