
...3TECH

TURBINE

Ti-940L

fr : Guide d'utilisation	03 - 21
en: User's guide	22 - 40
es : Manual de usuario.....	41 - 59
pt : Guia do utilizador	60 - 79
it : Manuale d'uso	80 - 98
de: Benutzerhandbuch	99 - 117
nl : Bedienungsanleitung	118- 137

TURBINE Ti-940L

fr

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TURBINE Ti-940L

en

Please carefully read this Instruction for Use before use.

TURBINE Ti-940L is manufactured from the best materials, with a meticulous finish, a pleasant design and reliable performance. The products are manufactured by combining German Technology with bearings from well-known brand bearing to create high performance. TURBINE Ti-940L can be compatible with Kavo® MULTiflex® LUX, Sirona® LED, W&H® Roto Quick Coupling, Bien-Air® Unifix® LED Coupling,, NSK® LED, Coupling, and will optimally meets the dentists' requirements.

1. Proper use

The medical device is only intended for dental treatment: removal of carious, cavities and crown preparation, removal of fillings materials, processing of tooth, restoration of surface and surgical procedure through the rotary processing by the compressed air to drive the milling cutter at high speed.

CAUTION

- The medical device is intended for use by a qualified knowledgeable user in dental workstation and laboratories.
- The medical device is intended for dental treatment. Any other type of use or alteration to the medical device is impermissible.
- Read carefully the Instructions for Use prior to using the medical device as it is always working at high speed during operation.
- During operation it is responsible for the user to protect the patient from danger and prevent contamination from the medical device.
- When beginning to run the medical device the user first should test run away far from the patient to confirm if it runs without any vibration, overheat and noise, otherwise the use must stop running.
- When beginning to run the medical device the user should take protective measures.
- During operation the user must stop running if an abnormality in the medical device occur.
- A damaged cartridge will make disturbing noise to cause hearing loss of the patient; it is impermissible to use the medical device with the

damaged cartridge.

- It is impermissible to disassemble or remake the medical device without authorization.
- It is always to prevent the medical device falling down the ground or being hit by something hard.
- When not in use a test bar (delivered with the medical device) should be inserted in the chuck system.
- As the medical device is non-sterile it should be sterilized in autoclaves before being used.

2. Product description

Type KAVO® Multiflex® LUX

Type Bien Air® Unifix® LED

Type W&H® Roto Quick


Type Sirona® LED




Type NSK® LED

3. Structure

- ① Push button
- ② Fiber optic
- ③ Grip profile
- ④ The turbine seat

- ⑤ BCooling spray nozzle
- ⑥ Cartridge
- ⑦ Nozzles needle
- ⑧ LED coupling

4. Features

TURBINE Ti-940L is manufactured from the best materials, with a meticulous finish, a pleasant design and reliable performance. The products find wide applications in dental treatment.

The cartridge design in structure according to international standard is assembled with the well-known brand bearings and the spindle system with run-out less than 0.008 and is accurately balanced for reliable rotation at high speed.

The patented TURBINE Ti-940L chucking system (Patent No.ZL2013 20057567.3) maintains more than 30N retention and a bur can be firmly chucked with accurate positioning.

TURBINE Ti-940L features quiet and vibration-free operation and offer a long service life and outstanding cutting efficiency.

TURBINE Ti-940L is designed with alternative to solid titanium body, brass body and stainless steel body to meet better user's requirements. The solid titanium body with non-allergic and biocompatible properties and in weight -- 30% less, more durable than the stainless steel or brass, and water-drop pattern on the body is to further improve grippable surface for maximum security in operation.

TURBINE Ti-940L offers head design with a 90° head angle or a 105° head angle for the user's alternative. Nozzle spray mode with 3-port or 4-port with requiring water amount more than 50 ml/min will offer effective cooling to the entire operation field.

TURBINE Ti-940L can offer the patented drive air pressure regulating system (Patent No.ZL2013 20060022.8) to maintain a stable supply of the pressure to the turbine to reach consistent rotation and hence better cutting effect.

Coupling system available in the Ti-940L turbines:

Type	Compatible with other brands speed coupling	Notice
KaVo® LED	KaVo® MULTIflex® LUX	K2 with spray ring
Sirona®	Sirona® LED coupler	With spray ring
W&H®	W & H® Roto Quick	With spray ring
Bien Air®	Bien-Air® Unifix® LED Coupler	B2 with spray ring
NSK®	NSK® PTL-CL coupler	With spray ring

5. Specifications and type

Head type	Torque head
Speed range	320,000 rpm - 350,000 rpm
Power	≥ 22 W
Retention power	≥ 30 N
Nozzle spray mode	4-port spray

Bur removal	Push button
Bur specifications	As per ISO1797-1
Coupling type	LED coupling
Voltage of LED Bulb	3 V
Illumination	≥ 25 000 LUX
Drive air pressure	250 ± 10% kPa (2.5 ± 10% kgf/cm ²)
Air consumption	40 NI/min ± 10%
Recommended water amount	> 50 ml/min

6. Assembly and removing the medical device

6.1 When LED coupling KAVO® Multiflex® is assembled

6.1.1 Assemble the LED coupling and attach the medical device

- Screw the LED coupling onto the turbine hose (Fig.6-1①).
- Tighten the screw cap of the turbine hose with the wrench (delivered with the medical device) (Fig.6-2②③).
- Assemble the medical device accurately on the LED coupling until it audibly locks in the medical device (Fig.6-2④).

6.1.2 Removing the medical device

Hold the LED coupling tightly; pull the medical device off while twisting slightly (Fig.6-3⑤).


Fig. 6-1


Fig. 6-2


Fig. 6-3

⚠ NOTE

Rotate the spray ring on the LED coupling in order to regulate the fraction of water when the LED coupling KaVo® Multiflex® is assembled. Recommended water amount should be more than 50 ml/min.

6.2 When LED coupling SIRONA® is assembled

6.2.1 Assemble the LED coupling and attach the medical device

- Screw the LED couplings onto the turbine hose (Fig.6-4①)
- Tighten the screw cap of the turbine hose with the wrench (delivered with the medical device) (Fig.6-4②③).
- Assemble the medical device accurately on the LED coupling until the coupling audibly locks in the medical device (Fig.6-5④).

6.2.2 Removing the medical device

Hold the LED coupling tightly; pull the medical device off while twisting slightly (Fig.6-6⑤).


Fig. 6-4


Fig. 6-5


Fig. 6-6

⚠ NOTE

Rotate the spray ring on the LED coupling in order to regulate the fraction of water when the LED coupling SIRONA® is assembled. Recommended water amount should be more than 50 ml/min.

6.3 When LED coupling W&H® Roto Quick is assembled

6.3.1 Assemble the LED coupling and attach the medical device

- Screw the LED coupling onto the turbine hose (Fig.6-7①).
- Tighten the screw cap of the turbine hose hang-tightly (Fig.6-7②).
- Assemble the medical device accurately on the LED coupling while sliding the lock ring (Fig.6-7③④).

6.3.2 Removing the medical device

- Hold the LED coupling tightly, slide the lock ring (Fig.6-8⑤).
- Pull the medical device off (Fig.6-8⑥).


Fig. 6-7


Fig. 6-8

⚠ NOTE

Rotate the spray ring on the LED coupling in order to regulate the fraction of water when the LED coupling W&H® Roto Quick is mounted. Recommended water amount should be more than 50 ml/min.

6.4 When LED coupling BIEN AIR® Unifix® is assembled

6.4.1 Assemble the LED coupling and attach the medical device

- Screw the LED couplings onto the turbine hose (Fig.6-9①).
- Tighten the screw cap of the turbine hose hang-tightly (Fig.6-9②).

- Assemble the medical device accurately on the LED coupling until the coupling audibly locks in the medical device (Fig.6-9③)

6.4.2 Removing the medical device

- Hold the LED coupling tightly, press two buttons on the side of the LED coupling (Fig.6-10④).
- Pull the medical device off (Fig.6-10⑤).


Fig. 6-9


Fig. 6-10

⚠ NOTE

Rotate the spray ring on the LED coupling in order to regulate the fraction of water when the LED coupling BIEN AIR® Unifix® is assembled. Recommended water amount should be more than 50 ml/min.

6.5 When LED coupling NSK® is assembled

6.5.1 Assembled the LED coupling and attach the medical device

- Screw the LED couplings onto the turbine hose (Fig.6-11⑩).

- Tighten the screw cap of the turbine hose hang-tightly (Fig.6-11②).
- Assemble the medical device accurately on the LED coupling until the coupling audibly locks in the medical device (Fig.6-11③).

6.5.2 Removing the medical device

- Hold the LED coupling tightly, slide the lock ring (Fig.6-12④).
- Pull the medical device off (Fig. 6-12⑤).


Fig. 6-11


Fig. 6-12

⚠ NOTE

Rotate the spray ring on the LED coupling in order to regulate the fraction of water when the LED coupling NSK® is mounted. Recommended water amount should be more than 50 ml/min.

7. Checking the pressure

Insert a test nanometer between the LED coupling and the medical device, check the following pressure:

- Drive air pressure (Standard head): $220 \pm 10\%$ kPa (2.2 ± 10 kgf/cm²)
- Drive air pressure (Torque head): $250 \pm 10\%$ kPa (2.5 ± 10 kgf/cm²)
- Return air pressure: ≤ 15 kPa (0.15 ± 10 kgf/cm²)
- Spray air pressure: $220 \pm 10\%$ kPa (2.2 ± 10 kgf/cm²)
- Spray water pressure: $200 \pm 10\%$ kPa (2.0 ± 10 gf/cm²)

NOTE

Make sure that the compressed air is dry, clean and uncontaminated in accordance to the Instructions for Use and set the pressure as per recommended value.

CAUTION

- If the compressed air is dirty and moist, it will cause premature deterioration of the rotor and premature wear of the bearings in the medical device.
- It is necessary to drain off water from the compressed air every day in order to protect the compressed air from moisture.
- Ensure a regular care and maintenance of the compressor and air-water filtering system.
- Ensure that the medical device is with dry, clean and uncontaminated air in accordance to ISO 7494-2.


Fig. 7

8. Inserting and removing the carbide bur or diamond bur

8.1 Insert the carbide bur or diamond bur

Carefully press the push button with your thumb and simultaneously insert the carbide bur or diamond bur (Fig.8-1①②).

8.2 Remove the carbide bur or diamond bur

Press the push button with your thumb and simultaneously pull out the carbide bur or diamond bur (Fig.8-1①③).


Fig. 8-1

⚠ NOTE

Only use the carbide bur or diamond bur that correspond to ISO1797-1 Type 3 and meet the following criteria:

- Shaft diameter: 1.59 – 1.60 mm
- Overall length max. 25 mm
- Shaft chucking length: at least 11 mm
- Blade diameter: max. 2 mm

⚠ CAUTION

When inserting the carbide bur or diamond bur

- Observe the Instructions for Use of the carbide bur or diamond bur and use them properly, especially rotational speed of the carbide bur or diamond bur should reach the required speed of the medical device, a lower speed of the carbide bur or diamond bur than required speed of the medical device is not permitted.
- Never use the carbide bur or diamond bur with worn, damaged or distorted shaft; otherwise they fall out during treatment to injure the patient. A

distorted shaft will reduce the bearings and rotor service life.

- Never use the worn and distorted blade of the carbide bur or diamond bur; otherwise they could cause premature deterioration of the rotor and premature wear of the bearings.
- Keep the carbide bur or diamond bur clean. A dirty or unclean shaft being inserted into the chucking system will cause unfirmed clamping.
- Insert the carbide bur or diamond bur to the bottom of the chucking system and pull on it to check that it is securely held and that the chucking system is okay.
- It is not permitted to use a short shaft of the carbide bur or diamond bur in the medical device with standard or torque head.
- Wear gloves or fingerstall to prevent an injury or infection.

When removing the carbide bur or diamond bur

- Do not touch rotating carbide bur or diamond bur.
- Never press push button while the carbide bur or diamond bur is rotating.
- Wear gloves or fingerstall to prevent an injury or infection.

9. Care and maintenance

9.1 Manual external cleaning and disinfection

- Brush off under flowing tap water with medium hard toothbrush
- Let the head upward and wash the body under flowing tap water to remove residual cement and composite or blood.
- Wipe down the medical device with a lint-free cloth immersed in cleanser
- Spray the disinfectant on a lint-free cloth, and thoroughly wipe down the medical device.

9.2 Manual internal cleaning and lubricating

- Use lubricant that is specified for the medical device to clean and lubricate the interior of the medical device: put the nozzle insert into the tube of the spray can and then insert the nozzle into the drive air tube of the medical device (Fig. 9-1), press the spray button for 2 – 3 seconds until the lubricant comes out from the head.
- Wipe off any leaking lubricant from the head with a lint-free cloth.

9.3 Chucking system cleaning and lubricating

- Use lubricant that is specified for the medical device to clean and lubricate the chucking system: put the nozzle insert into the tube of the spray can and then insert the nozzle to the chucking System (Fig. 9-2), press the spray button for 2 – 3 seconds until the lubricant comes out from the head.
- Wipe off any leaking lubricant from the head with a lint-free cloth.


Fig. 9-1


Fig. 9-2

9.4 Manual drying

Blow off the exterior and interior with the compressed air until no water drops are visible

9.5 Sterilization

- Place the medical device individually a in the sterilization bag, and put into a steam sterilizer (autoclave)
- Sterilize in autoclave with by performing 3 pre-vacuum cycles for at least 18 minutes at 134°C (273°F).

NOTE

- Improper servicing and care of the medical device will cause the premature wear and malfunction, thus will reduce the medical device life, it is requested to service and care the medical device according to ISO 17664 to proceeding with the procedure: cleaning, disinfecting, lubricating and sterilizing.
- It is not permitted to wash the head of the medical device.
- In the procedure of cleaning the exterior of the medical device it is permitted to spray the disinfectant on a lint-free cloth, then thoroughly wipe down the medical device.
- In the procedure of cleaning the interior of the medical device if any dirty leaking lubricant come out, it is requested to repeat the lubricating procedure.
- If it is necessary attach the medical device to the LED coupling to the compressed air hose to run for a 10-20 seconds to distribute the lubricant to the cartridge of the medical device more effectively
- Hold the medical device firmly to prevent the medical device falling down from your hand due to excessive spray pressure.
- In the procedure of chucking system cleaning and lubricating hold the spray can upright.
- Carry out chucking system cleaning and lubricating at least once a week.
- Always hold the spray can upright, not upside-down

CAUTION

- In cleaning procedure the cleanser should be neutralization agent, never use the cleanser or disinfectant containing chlorine.
- Do not immerse the head of the medical device or the medical device to a disinfectant bath, and it is not permitted to clean the medical device in the ultrasonic unit.
- Autoclaving can be operated in a steam sterilizer that corresponds with EN13060 B/ISO17665-1 Class B or S and is specified for the dental turbines.
- Autoclaving maximum temperature does not exceed 138°(280°F)
- Remove the medical device from the steam sterilizer and store in a dry, dark and cool space protected from dust with minimum exposure to germs.
- Do not immerse the medical device in cold water to speed up the cooling process after removing from the steam sterilizer.

10. Storage

- Reprocessed medical device should be stored in a dry, dark and cool space protected from dust with minimum exposure to germs, corrosive gas and far away from the heat source.
- Comply with the expiry date of the sterilized items.

11. Troubleshooting

Failure	Root cause	Troubleshooting
Excessive noise	Worn bearings	Replace new bearings
Lower rotating speed	Broken bearings Insufficient drive air pressure	Replace new bearings Re-set air pressure
Insufficient cutting power	Broken bearings Insufficient drive air pressure	Replace new bearings Check and set air pressure
Non-rotation of rotor	Broken bearings	Replace new bearings
No spraying Insufficient spraying	Blocked water passage or Insufficient water and air pressure	Clean water channel with the nozzle needle Check and set water and air pressure
Less retention of bur	Non-standard bur Dirt in chucking system	Use the standard bur Clean and lubricate chucking system
Vibration of bur rotation	Non-standard bur Broken bearings	Use the standard bur Replace new bearings
No lighting of fiber optic	Broken LED bulb	Replace news LED bulb

12. Warranty terms and conditions

Our company provides the user with a 12-month warranty on all our product range, with the exception of ball bearings, from the date of invoice. Maintenance during the warranty period is the responsibility of the customers.

Our company is not liable for damages or injuries resulting from:

- Excessive use.
- Any improper handling of the product or any modification made to the product by persons not authorized by our company.
- Failure to comply with the instructions for installation, use and maintenance of the handpiece.
- Any chemical, electrical or electrolytic damage caused by improper autoclaving and/or storage.
- Inappropriate working pressure.

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Shanghai International Holding Corp. GmbH
Tealth Foshan Medical Equipment Co., Ltd