**Shanghai International Holding Corp. GmbH**  
Eiestrasse 80, 20537 Hamburg, GERMANY

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## SANDBLASTING MACHINE WITH WATER COOLING SYSTEM

Read the instructions carefully before using them to improve working efficiency and reduce failure.

### 1. Principle

Sand powder particles are sprayed through high-pressure airflow (18 m/s) are used to roughen surfaces or to remove residues of various materials. MICRO SANDER AW+ uses a special medium for aluminum oxide sandblasting powder or sodium bicarbonate powder (containing Glycine, the same below). MICRO SANDER AW+ adopts a built-in water supply pipe and deep embedded tungsten steel nozzle, which is compact, light, and easy to grip. It is an ergonomic design. During the use process, the water cools down the powder and effectively inhibits powder spatter.

### 2. Product application

1. Use of alumina sandblasting powder for surface cleaning and micro roughening before bonding (Fig. 1) such as hole preparation / minimally invasive treatment, removal of sundry/coarse bonding face/removal crown in the groove, residual cement on the bridge and metal casting repair body defects detection.
2. The use of sodium bicarbonate micro powder can be applied to a variety of clinical applications (Fig. 2) such as removing staining and plaque, tobacco stains and biofilms to reduce the risk of dental caries. At the same time, after using ultrasonic cleaning, sandblasting can also play a role in polishing and repairing the surface of the teeth.
3. When applied in the mouth, given the assistance of water sources during sandblasting, the pressure directly acting on the teeth can be reduced, the cleanliness is efficient, the target is accurate, and the patient does not feel dry due to dust interference.

en

Fig. 1

Fig. 1 shows the situation observed when the prepared surface of the hole is enlarged 400 times after using the aluminum oxide powder microroughness, and it can be clearly seen that the entire surface is evenly rough.

Fig. 2 is a comparison of the use of sodium bicarbonate powder before and after oral cleansing. Comfortable, efficient, safe, without damage to the dentin, enamel and orthodontic support groove, protects the teeth and permanent teeth.


Fig. 2

### 3. Product Installation

1. The connector of the MICRO SANDER AW+ is directly connected to the dental unit handpiece tube (Fig. 3). (Including turbines and other water supply and gas equipment, the same below).
2. Turn off the water and air supply to the dental unit, remove the MICRO SANDER AW+ powder bottle, pour it into the sand powder, where the sand powder must not exceed the maximum capacity indicated by the bottle body, and then screw the bottle into the handle body and check whether it is tightened (Fig. 4).
3. The nozzle is aligned with the spittoon, the supply switch on the dental unit is first turned on, the foot switch is stepped down, and the starting key on the handle of the MICRO SANDER AW+ is lightly pressed. After the sand powder is ejected and adjusted evenly, then turn on the water switch on the dental unit and adjust the amount of water. Try to spray 3 seconds first, confirm that the sand powder injection achieves the desired effect, and then use the target (Fig. 5). (Water and air adjustments are done on the dental unit. This process does not need to start the atomization switch).

Fig. 3

Fig. 4


Fig. 5

#### Warm tips:

- Use of alumina sandblasting powder: The standard pressure on the dental unit should be adjusted to 3.5-5 .5 kgf/cm. It is recommended to use < 50  $\mu\text{m}$  high-purity alumina sandblasting powder and not to block for moisture. The operator shall wear a mask and goggles and operate in a relatively sealed transparent dust collection box. Do not spray the nozzle on a non-sprayed target to avoid loss or danger. It is recommended to use rubber dam and strong suction system to assist in sand blasting.
- Use of sodium bicarbonate micro powder: The standard pressure on the dental unit should be adjusted to 2.5-4 .5 kgf/cm. It is recommended to use < 65  $\mu\text{m}$  sodium bicarbonate or Glycine micro powder and not to block for moisture. During operation, the nozzle should be 3-5 mm away from the surface of the tooth, and a small circle motion jet should be made at an angle of 30-80 degrees. It is strictly forbidden to spray sand powder directly into the soft tissue of the mouth, gums, or to use it in patients who are not suitable for cleaning the tooth surface with sodium bicarbonate powder. After the treatment is completed, adjust the amount of water on the dental unit or turbine to the maximum surface of the rinsing tooth. (The operation should be performed by a dentist and a dental health care worker. Persons who have not obtained the relevant qualifications are not allowed to use it).
- **At the end of the operation, the MICRO SANDER AW+ should be removed to a non-sandblasting target. The water supply should be shut down first, and then the gas supply should be closed so that the accumulated water from the water hole will not seep into the sand nozzle, causing the sand nozzle to block. The remaining sand powder shall not be reused, and the empty powder bottle shall be tightened to the sand-cleaning machine handle body, and the empty spray shall be carried out when the water supply is closed until the handle pipe is emptied without residue.**

### 4. Operation Tips

1. If the MICRO SANDER AW+ accidentally falls, check whether the powder bottle has cracks and whether the nozzle is blocked before re-use. If there is damage, it should be replaced immediately. The water supply and air switch must be closed when not used for a short period of time.
2. The amount of sand powder shall not be higher than the maximum range defined by the powder bottle. When the capacity of the powder bottle is less than 1/5, the water supply switch shall be closed immediately to make up the sand powder in time.
3. If the target is sprayed with alumina sand powder in the dust collection box, the water supply should be closed at this time. The target should only be dry sprayed to avoid the danger of the dust collection box being exposed to water leakage.
4. The metal part of the MICRO SANDER AW+ can be sterilized at 134°C. Other detachable parts such as powder bottles and O-rings and related gaskets are recommended for alcohol wiping disinfection.

**REMARK: DO NOT PROCESS THE SAND-SPRAY CLEANING MACHINE AS A WHOLE FOR HIGH TEMPERATURE STERILIZATION, IT SHOULD BE DISASSEMBLED AND SELECTIVE MAINTENANCE DISINFECTION APPLIED.**

## 5. Troubleshooting

### If blockage occurs:

1. Remove the powder bottle, turn off the water supply, plug the nozzle out of the sand hole with a soft adhesive pad, set the air pressure to maximum (6.0 kgf/cm), step on the pedal switch, and press the start key for a few seconds. Release the nozzle out of the sand hole, use the soft adhesive pad to alternately block the air hole and sand hole in the powder bottle, and drain the nozzle clean.
2. Remove the nozzle with the attached wrench and reverse the nozzle to shake off the powder inside. Before loading the nozzle into the handle, please step on the foot switch first, gently press the start key to spray a little of the residual sand in the handle until the handle is fully dredged, and then the nozzle is twisted into the handle.
3. If none of the above methods can be solved, please contact your distributor.

## 6. Guarantee

The MICRO SANDER AW+ is guaranteed for 1 year. Attachment and human damage are not covered by the free warranty. In the event of unauthorized self-disassembly the manufacturer reserves the right to no longer provide free warranty. Unclogging not included in the warranty. The manufacturer reserves the right to carry out product improvement without prior notice.

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