Shanghai International Holding Corp. GmbH
Eiffestrasse 80, 20537 Hamburg, GERMANY

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Date de fabrication

LED SURGICAL HANDPIECE WITH EXTERNAL IRRIGATION

CAUTIONS

Please read this instruction before use and file for future reference.

Prior to each use, all instruments and accessories must be inspected for proper operation.

When using this product, pay attention to protect the patient at all time.

This product should not be rebuilt. Replacing bearings only could result in unexpected failure and patient injury due to different factors of new vs. old components.

End of product life, the metal can be handled by separation and recycling disposal.

Products storage and transportation should be carried out in environment with the following features: dry, stable (no collision), clean, normal temperature, normal pressure, away from acid gases and other harmful chemical substances, and humidity should be less than 90%.

Maintenance should be made by professionals.

This product should be used with caution for the patients with bleeder disease, pacemaker, and pregnant woman or children.

en

1. Product description

This medical device is:

- intended for dental treatment only. All other types of use or alterations to the product are not permitted and can be hazardous. The medical product is designed to be used with the corresponding heads in the following applications: surgery e.g. bone augmentation, sinus lift, dental extraction, implantology and maxillo-facial surgery.

- A medical device according to relevant national statutory regulations. This product can be connected to motors with a connection in accordance with EN ISO 3964. During the operation, Sterile physiological coolant is required, and water spray on the unit must be turned off.


Fig. 1

2. Technical Specification

Drive speed	Max. 40,000 rpm
Bur application	ISO 1797 ϕ 2.334-2.35 mm
Speed ratio	1:1
Connection	ISO 3964
Noise	\leq 70 dB


Fig. 2

3. Operation

3.1 Attach the handpiece

Place the medical device on the motor clutch and lock it into place. The latch must lock into place audibly (Fig.3).
Make sure that the handpiece is securely affixed to the motor.


Fig. 3

3.2 Remove the handpiece

Unlock the medical device from the motor coupling by twisting it slightly and then pulling it along its axis (Fig. 4).


Fig. 4

CAUTION

When motor is running, do not remove the handpiece.
Stop operation if the handpiece is getting warm.

3.3 Insert the handpiece drill bur

- Rotate the tension ring in the direction of the arrow to the backstop and insert the surgical bur into the chuck.
- Turn the tension ring back into its initial position.
- Check that the surgical bur is securely attached by pulling on it.


Fig.5

CAUTION

- Do not use unauthorized surgical burs.
- Make sure that the surgical bur is inserted till the end.
- Do not run the motor when the tension ring is not back into its initial position.
- Keep a bur inserted in the handpiece even when it is not used.

4. Cleaning

- Clean the handpiece after each use.
- Clean the surface with a toothbrush and wet wipes.
- Clean the spray tube with nozzle needle to free the water passage.

⚠ CAUTION

- Do not soak this product into chlorine disinfection.
- Do not clean this product with ultrasonic cleaner.

5. Maintenance

Lubricate the handpiece at least one time each day, 1-2 seconds each time.


Fig.6

6. Sterilization

Sterilize this product in a steam sterilizer in accordance with ISO 17665-1, 134°C, for at least 18 minutes.

7. Terms and conditions of warranty

The manufacturer provides the end customer with a warranty of proper function and guarantees zero defects in respect of material and processing for a period of 12 months from the date of the invoice, subject to the following conditions:

In case of justified complaints, the manufacturer will honor its warranty with a free replacement or repair. Other claims of any kind whatsoever, in particular with respect to compensation, are excluded. In the event of default, gross negligence or intent, this shall only apply in the absence of mandatory legal regulations to the contrary.

The manufacturer shall not be liable for defects and their consequences that have arisen or may arise from natural wear, improper handling, cleaning or maintenance, noncompliance with operating, maintenance or connection instructions, calcination or corrosion, contaminated air or water supplies or chemical or electrical factors deemed abnormal or impermissible in accordance with instructions for use. The warranty granted does not usually extend to lamps, optical fibers made of glass and glass fibers, glassware, rubber parts, and the colorfastness of plastic parts.

All liability is excluded if defects or their consequences originate from manipulations or changes to the product made by the customer or a third party that is not authorized by the manufacturer.

Warranty claims will only be accepted if the product is submitted along with proof of purchase in the form of a copy of the invoice or note of delivery. The dealer, purchase date, type, and serial number must be clearly evident from this document.

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